

Agilent Seahorse XFp Analyzer
Specifications

Item	
Number of Assay Wells	6
Microchamber Volume	2 μ L
Controller	Integrated full-color touch-screen interface to design basic protocols and run imported templates.
Software	Create and run standard assays on controller OR design on Wave desktop, transfer to Controller, return to desktop for analysis.
Best for	<ul style="list-style-type: none">• Pairwise comparisons• Phenotyping single samples• Precious biomaterial• Patient-derived samples
Key Advantages	<ul style="list-style-type: none">• Easy to set up and run via streamlined interface• Standardized protocols for XF Assays• Validated for assay temperatures 16-40 °C (8-20 °C lower ambient temp required)
Plate Format	8-well microplate