

Pacific Biosciences PacBio Sequel System
Specifications

Item	
Movie Time	30 mins - 10 hours (Maximum Data Acquisition time: 72 hours)
Loading Method	<ul style="list-style-type: none"><li data-bbox="579 353 1098 387">• Diffusion Loading (for library <7500 bp)<li data-bbox="579 398 1098 432">• MagBead Binding (for library >7500 bp)
Pre-Extension Time for Amplicons	30 - 120 mins (With the use of Pre-Extension Feature, longer average polymerase read length can be achieved.)