

Nikon Fluorescence Microscope (Eclipse Ti Motorized Inverted Microscope)
Specifications

Item	
Objectives	<ul style="list-style-type: none">• 4x/ 0.3• 10X/ 0.3• 20X/ 0.45• 40X/ 0.95 OIL• 60X/ 1.4 OIL• 100X/ 1.45 OIL
Filter Sets	<ul style="list-style-type: none">• DAPI• FITC• Rhodamine• Alexa 568• Cy5• Transmitted light
Camera	Andor Zyla 4.2 sCMOS, 4.2 Megapixel
Stage	Manual stage
Sample Insert	Universal
Software	Nikon NIS Element